

## IPS Contrast Spray Chairside Refill Cannulas

IPS Contrast Spray Chairside

JDI003



## IPS Ceramic Etching Gel Red

Used for the creation of retentive bonding surfaces for the composite cementation technique. It enhances the bonding effect between the composite cement and the ceramic bonding surface.

JDI004



## IPS Empress Veneer Kit

Contains the most popular Wash pastes and layering materials for the veneering of IPS Empress restorations (Esthetic & CAD).

JDI242



## IPS Neutralising Powder

Used for the creation of retentive bonding surfaces for the composite cementation technique.

JDI005



## IPS e.max ZirCAD LT Colouring Liquid Starter Kit

IPS e.max ZirCAD Colouring Liquids are ready-to-use, aqueous metallic salt solutions for colouring unsintered restorations made of IPS e.max ZirCAD MT BL and LT.

JDI807



## IPS e.max ZirCAD MT Colouring Liquid Starter Kit

IPS e.max ZirCAD Colouring Liquids are ready-to-use, aqueous metallic salt solutions for colouring unsintered restorations made of IPS e.max ZirCAD MT BL and LT.

JDI808

